


ARUM



iPen

Lab driver set

[Screw driver tips]



- Place the model scanbody and the abutment
- 10 types of tips (compatible with the various implant systems)



Ball Torx



Ball Hex



Hex 1.0



Hex 1.2



Hex 1.25



Hex 1.27



Hex 1.3



Torx




Star



Cross

[Analog insertion tips]



- Place the lab analog to 3D printed model
- Easy and safe to place the lab analog
- Prevent damage of the lab analog internal thread
- 5 types of tips (compatible with the various implant systems)



M1.4*P0.30



M1.6*P0.35



M1.8*P0.35




M2.0*P0.40



M2.5*P0.45

[Torque adapter]



- For clamping Attachment and Premill for connectorless milling (with 0.4Nm torque)
- Prevent screw loosening during milling



Screw driver tips line-up

Size	Colour	Compatible Implant system
Hex 1.0	Blue	DENTIUM® NR Line / DENTSPLY® Ankylos® C/X Implant etc.
Hex 1.2	Purple	ADIN® TOUAREG™ S&OS / ANTHOGRY ANTHOFIT / BIOMET 3i® Certain / BIOMET 3i® EXTERNAL HEX / BTI INTERNAL® / COWELLMEDI® INNO Internal / DENTIS® i-Clean / DENTSPLY® Xive® / DIO®UF / Global D® Tekka® In-Kone / MEGAGEN® AnyONE / MEGAGEN® Anyridge® / MEGAGEN® Anyridge® / MEGAGEN® Extra EZ POST / MEGAGEN® EZ Plus / MEGAGEN® Rescue / NEOBIOTECH® / OSSTEM® GS(TS) / OSSTEM® SS / OSSTEM® US / SHINHUNG® / SIC Invent®/SOUTHERN IMPLANTS MSc EXTERNAL HEX / WARANTEC® Oneplant etc
Hex 1.25	Black	BEGO® Semados / BIOHORIZONS® External® / BIOHORIZONS® Internal® / DENTIS® s-Clean / LASAK Bioniq® etc.
Hex 1.27	Red	ANTHOGRY AXIOM / Alpha-Bio Tec® Internal Hex / ASTRA TECH™ OsseoSpeed™ TX / CAMLO® CAMLO® / CAMLOG® CONELOG® / CORTEX™ / DEEP® / DENTIUM® SimpleLine / DENTIUM® SuperLine / EBI®/ IMPLANT DIRECT® Legacy™ / MIS® C1 / MIS® Internal Hexagon / ZIMMER® Spline A / ZIMMER® Spline B / ZIMMER® Tapered Screw-Vent® etc.
Hex 1.3	Green	ASTRA TECH™ OsseoSpeed™ EV / CLC Conic / DYNA Pushin Octalock® KENTEC® / PLATON® Type IV Pro / Sweden & Martina etc
Torx 1.7	Brown	STRAUMANN® Bone Level ® / STRAUMANN® ConicalBone Level ® / STRAUMANN® synOcta® etc.
Torx 1.5 (Star)	Yellow	ADIN® CLOSEFIT™ / NOBEL BIO CARE® Active™ / NOBEL BIO CARE® Branemark ® / NOBEL BIO CARE® Replace Select™ / NEOSS Pro Active® etc.
Cross	Grey	THOMMEN SPI® etc
Ball Torx	Orange	1. Most of angled screws of ARUM Ti-base (screw head size should be more than 2.1mm) 2. ALL IMPLANT SYSTEM FOR DISK ANGULATED SCREW CHANNELS
Ball Hex	White	Some of angled screws of ARUM Ti-base (screw head size should be less than 2.1mm) *Please visit our website to see the list.

Analog insertion tips line-up

Size	Colour	Compatible Implant system
M2.0*P0.40	Cocoa	Osstem GS Regular / Dentium Super Line / Megagen Anyone / Dentis S-clean / Straumann Tissue Level / etc.
M1.8*P0.35	Light Green	Zimmer TSV / Megagen Anyridge / Biohorizon Internal / Bego Semados / etc.
M1.6*P0.35	Pink	Straumann Bone Level / Osstem GS Mini / Nobel Biocare Active Mini / Dentsply Xive / Neodent GM / etc.
M1.4*P0.30	Sky Blue	Straumann Bone Level 2.9 / Dio UF Narrow / Megagen Mini / Astra Yellow / Most of Multi-unit / etc.
M2.5*P0.45	Ivory	Osstem US Wide / Nobel Biocare Branemark WP / etc.

Torque adapter

Torque	0.4 Nm
--------	--------

** Aware of the colour of analog tip

** Contact us for a complete list of implant compatibility

